



Microplate reader operating procedure

Shuning Guo¹, Zhujun Wei¹

¹2021 iDEC NEFU_China

Oct 05, 2021

1 *Works for me*

Share

dx.doi.org/10.17504/protocols.io.bysspwee

2021 iDEC NEFU_China



Shuning Guo

ABSTRACT

Use Tecan Spark® multimode microplate reader to measure absorbance or fluorescence intensity of green fluorescence protein.

SPARK

Microwell plate reader

TECAN

SPARK



DOI

dx.doi.org/10.17504/protocols.io.bysspwee

PROTOCOL CITATION

Shuning Guo, Zhujun Wei 2021. Microplate reader operating procedure.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.bysspwee>



KEYWORDS

multimode microplate reader

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED


Oct 05, 2021

LAST MODIFIED

Oct 05, 2021

MATERIALS TEXT

Tecan Spark® multimode microplate reader



96-well plate

SAFETY WARNINGS

Please read the instructions carefully before use.

Please turn off the power after use.

BEFORE STARTING

1. Click on "instrument" on the left, click on "temperature", and check "temperature control" to control the real-time temperature.
2. In the "Action" bar of the left column, you can drag it to the process list under the 96-well plate, drag up and down to sort, click the trash can to delete the option. You can design your own processes that meet your needs.

- 1 Turn on the power and warm up for 10 minutes.
- 2 Open the program and select "New" to create a new process.
- 3 Select "Absorbance" in the left column. There is a schematic of a 96-well plate in the middle of the screen.
- 4 You can select the hole to be used according to the spotting hole, and the selected hole will turn blue.
- 5 Put in the 96-well plate, set the parameters, click "start" to measure the absorbance, and then remove the 96-well plate.

- 6 Wait a moment, it will generate an Excel file with results and various parameters. Save and copy the file to the u disk and process the data.