



Nucleic acid & protein electrophoresis V.2

Shuning Guo¹

¹2021 iDEC NEFU_China

Version 2 ▼


1 Works for me

Share

Oct 07, 2021

dx.doi.org/10.17504/protocols.io.byu4pwyw

2021 iDEC NEFU_China

 Shuning Guo

ABSTRACT

This protocol concludes two types of the electrophoresis used to detect target DNA or protein.

DOI

dx.doi.org/10.17504/protocols.io.byu4pwyw

PROTOCOL CITATION

Shuning Guo 2021. Nucleic acid & protein electrophoresis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.byu4pwyw>
Version created by Shuning Guo



KEYWORDS

electrophoresis, SDS-PAGE

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 07, 2021

LAST MODIFIED

Oct 07, 2021

PROTOCOL INTEGER ID

53884

MATERIALS TEXT

Agarose electrophoresis:

Microwave

Electrophoresis apparatus

Electrophoresis pool

Gel bed

DdH₂O

Agarose

50X TAE

Gold view

6 X orange loading buffer

SDS-PAGE:

Microwave

Centrifuge

Electrophoresis apparatus

Electrophoresis pool

Glass plate

casting stand

ddH₂O

30% Acr- Bis(29:1)

1.5 mol/L Tris-Gly (pH 8.8)

10% SDS

10% Ammonium persulfate

TEMED

1.0 mol/L Tris-Gly (pH 6.8)

5× Tris-Gly electrophoresis solution (Tris-Base 15.1g, Glycine 94g, SDS 5g, pH=8.3)

coomassie blue staining solution (400ml, coomassie bright blue r-250 0.4g, isopropyl alcohol 100ml, glacial acetic acid 40ml, ddW260ml, filtered)

Destaining solution (500ml, glacial acetic acid 50ml, anhydrous ethanol 75ml, distilled water 375ml).

SAFETY WARNINGS

It is necessary to wear latex gloves and mask when preparing the gel because most of the reagents are toxic.

All the things used in the next steps should be washed after use and can not be touched without wearing gloves.

Disposable plastic gloves are not allowed to use.

BEFORE STARTING

Prepare 50 x TAE, 1.5 mol/L Tris-Gly (pH 8.8), 1.0 mol/L Tris-Gly (pH 6.8), 5× Tris-Gly electrophoresis solution, coomassie blue staining solution, destaining solution before start.

- 1 Choose suitable electrophoresis method depends on the type of the sample.

Step 1 includes a Step case.

Agarose gel electrophoresis

SDS-PAGE

Agarose gel electrophoresis

DNA was detected by agarose gel electrophoresis.

- 2 Weigh appropriate agarose depends on the concentration of the gel (1% agarose gel for detection and 1% or 2% for gel extraction).
- 3 Add 1X TAE to a conical flask. Need to prepare 1X TAE with 50X TAE.

Extra 2~3ml of 1X TAE was strongly recommended to be added to avoid the reduction of solution during heat.

- 4 Heat up by microwave until the solution is homogeneous.
- 5 Cool at room temperature for 3~5 min.
- 6 Add 1 μ l Gold view into solution and mix well when the solution is about 60°C.
- 7 Put the solution into bed for polymerize, make sure "comb" is well placed and the solution is balanced.
- 8 Wait about 20 min to let the gel completely concretes.
- 9 Mix the sample with 6 X orange loading buffer and load the sample into the sample holes.

Remember to load the DNA marker to the sample hole.

- 10 Put the bed with gel into the electrophoresis chamber.

11 Set the voltage of electrophoresis (80V~150V) and begin to run.

12 Stop running when the front indicator reach about 3/4 length of the gel.

13 Use

Gel Doc XR+ Gel Documentation System
Gel Documentation System

Bio-rad Laboratories 1708195 [🔗](#)

to observe the gel.