



Overlap & Gibson ligation

Shuning Guo¹

¹2021 iDEC NEFU_China

Oct 05, 2021

1 *Works for me*

Share

dx.doi.org/10.17504/protocols.io.byrspv6e

2021 iDEC NEFU_China



Shuning Guo

ABSTRACT

This protocol is used to ligate two pieces of DNA together without digesting the fragment with restriction endonucleases.

DOI

dx.doi.org/10.17504/protocols.io.byrspv6e

PROTOCOL CITATION

Shuning Guo 2021. Overlap & Gibson ligation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.byrspv6e>



KEYWORDS

Overlap, Gibson assembly

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 04, 2021

LAST MODIFIED

Oct 05, 2021

PROTOCOL INTEGER ID

53778

MATERIALS TEXT

DNA fragments

Primers

ClonExpress II One Step Cloning Kit (Vazyme)

2×High Fidelity Master Mix (MCLAB)

Nanodrop

Thermocycler

Water bath

DdH₂O

SAFETY WARNINGS

Please wear gloves for the experiment, don't try to touch the lid after PCR program initiation.

BEFORE STARTING

Set up a small box with ice, put DNA and enzymes on it.

Prepare the water bath to 37°C to have Gibson assembly.

- 1 Ligation of two DNA fragments by using cases below.
Step 1 includes a Step case.

Overlap PCR

Gibson Assembly

Preparation of linearized vectors

step case

Overlap PCR

- 2 Select an appropriate cloning site on the vector that will be linearized.
- 3 Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes or by reverse PCR.

PCR of the inserts DNA fragments

- 4 Amplify the insert DNA fragments with homologous sequences (for homologous recombination) of vector-upstream or -downstream by PCR using high fidelity DNA polymerase.

Calculate amount and ratio of linearized vectors and Inserts

- 5 Detect DNA concentration of linearized vectors and inserts by Nanodrop.
- 6 Calculation of the amount of vectors:
Molar ratio of vector to insertion is 1:1

Recombination & PCR

7 Set up the following reaction on ice (50 μ l):

A	B
Forward Primer (10 μ M)	1 μ l
Reverse Primer (10 μ M)	1 μ l
Fragment1(vector)	X
Fragment2(insertion)	Y
2 \times High Fidelity Master Mix (MCLAB)	25 μ l
ddH ₂ O	Add to 50 μ l

The primer is used to amplify recombinant DNA fragment/circular DNA.

8 Program the thermocycler as follows:

A	B
Temperature	Time
95/98 $^{\circ}$ C	5 min
95/98 $^{\circ}$ C	30 s
T _m -3~5 $^{\circ}$ C	30 s
72 $^{\circ}$ C	1kb/min
72 $^{\circ}$ C	5~10 min
16 $^{\circ}$ C	∞

Repeat 30 times in 3-5 steps

9 Use the palm centrifuge to mix the solution in PCR tube.

10 Put the PCR tube into the thermocycler and Run the program.

11 Using agarose gel electrophoresis to confirm if correct construct was present.